



# Accelerating marker-assisted selection - from dirt to data

Dr. Jonathan Curry – Senior Scientist - Genomics

A photograph of a wheat field under a bright blue sky with scattered white clouds. The wheat stalks are golden-brown and in sharp focus in the foreground, with the rest of the field and sky blurred in the background.

Science  
for a safer world



# Introduction

- Who are LGC?
- Where do SNPs grow?
- How might they be harvested for breeding?
- You reap what you sow – applying technology

# LGC Genomics is a part of the LGC Group



- LGC established in 1842 as a scientific testing laboratory for the UK government
  - Privatised in 1996
- Retains its role as UK National Measurement Institute
  - Standard bearer for chemical and bio-analytical measurement
- Global presence
  - 33 locations
  - 2,000+ staff
  - >3-fold growth since 2004

## LGC Health Sciences

Leader in analytical testing services for food, drug, and consumer products



## LGC Genomics

Genomic solutions for DNA/RNA extraction, genotyping, and sequencing



## LGC Forensics

The world's largest supplier of products and services for forensic analysis



## LGC Standards

Production and supply of chemical and biochemical reference materials



## LGC Science & Technology

Home of the UK National Measurement Institute




**LGC Genomics  
Division - Hoddesdon,  
Herts, UK**

**LGC Genomics  
Division - Berlin,  
Germany**

**LGC Genomics  
Division - Beverly, MA  
USA**

**LGC locations**






# LGC Genomics Labs - Huddesdon






# SNP genotyping project process – Kraken™ integration check points



**Kraken™ integrated SNPlines checks at vital points of the process ensure the correct sample and project are always tied.**

# Distribution of Life in Kraken



■ Plant ■ Animal ■ Fish ■ Insect ■ Bacteria / Yeast / Protozoan ■ Human

# Where do SNPs grow?

A photograph of yellow flowers, likely rapeseed, against a clear blue sky. The image is framed by a teal border with a rounded corner on the right side.

Science  
for a safer world




# High-Throughput Marker Discovery

Discovery of valuable traits using Next-generation sequencing




- Using High-throughput sequencing to survey genomes in bulk
- Markers are used to identify differences between genomes.



- Use markers when cross breeding varieties.
- Once markers are selected they can be run over hundreds / thousands / millions of samples.
- Reduction of time to market – reduces cost / increases profitability



# KASP genotyping landscape



# Trait Identification




The reduction in cost of NextGen sequencing has meant a dramatic increase in possible loci to target

## Position of desired trait marker



## Fine mapping with KASP genotyping




## Design flexibility allows more markers to be targeted to saturate regions

By using more SNP markers (higher resolution mapping) you will find desired trait faster

KASP assays allow higher resolution, flexible mapping of organisms

# High Throughput Marker Assisted Selection

LGC



**How might they be harvested for breeding?**

A photograph of yellow flowers, likely rapeseed, against a clear blue sky. The image is framed by a rounded teal border on the right side.


**Science**  
for a safer world




# KASP Components




A) Assay mix




B) Master mix



Legend



C) DNA template (sample)



- Assay mix - contains two allele specific primers and common reverse primer(s). This is often referred to as a KASP Assay.
- Universal Master mix - contains fluorophores, quenchers and *Taq* polymerase.

# KASP – flexible for all platforms

KASP is very simple and scalable...

- Anyone who has a pipette, PCR machine and a way of reading FRET can run KASP!!
- System agnostic
- LGC Genomics service labs can run your project for you.
- Bi-allelic SNP data is general purpose marker currency – easy to exchange and work with.

# Flexible design



- KASP is PCR
  - Can be designed as such
  - Can be optimised as such (Temp / Betaine / DMSO/  $Mg^{2+}$ )
  - 90% assay conversion
- Choose loci with a fair chance of designing primers.
- i.e. not:
  - Non-unique and highly spread throughout genome (although this can be solved)
  - Long homo-polymer repeats i.e. AAAAAA
  - VNTRs / micro-satellites – although Hybeacons can be used here.
  - Copy number




# Solutions for homology and polyploidy



Want to save time and money on breeding / mapping projects with KASP?


- Answer?
  - BLAST, BLAST and BLAST some more!!!!!!
- Very simple to do with available databases for many common crops.
- Turn off complexity filters and run Mega-BLAST then BLASTN.
- Even partially sequenced genomes can give some information about choice of markers to design.

# Solutions for homology and polyploidy





# Polyploid data – allele dosage



## Heterozygotes.

In a tetraploid organism these samples could have any of the following genotypes:

1:1:1:2

1:1:2:2

1:2:2:2

- The central set of data points are from known heterozygotes.
- These are wheat crosses with clear allele dosage effect seen centrally.
- If using this kind of experiment then well normalised DNA is crucial.

**You reap what you sow – applying  
technology**

A photograph of yellow flowers, likely rapeseed, against a clear blue sky. The image is framed by a teal border with a rounded corner on the right side.

**Science**  
for a safer world






EMBL



TEXAS A&M  
UNIVERSITY



UNIVERSITÄT LEIPZIG



## Discovery and development of exome-based, co-dominant single nucleotide polymorphism markers in hexaploid wheat (*Triticum aestivum* L.)

Alexandra M. Allen<sup>1,\*</sup>, Gary L. A. Barker<sup>1</sup>, Paul Wilkinson<sup>1</sup>, Amanda Burridge<sup>1</sup>, Mark Winfield<sup>1</sup>, Jane Coghill<sup>1</sup>, Cristobal Uauy<sup>2</sup>, Simon Griffiths<sup>2</sup>, Peter Jack<sup>3</sup>, Simon Berry<sup>4</sup>, Peter Werner<sup>5</sup>, James P. E. Melichar<sup>6</sup>, Jane McDougall<sup>7</sup>, Rhian Gwilliam<sup>7</sup>, Phil Robinson<sup>7</sup> and Keith J. Edwards<sup>1</sup>

<sup>1</sup>*School of Biological Sciences, University of Bristol, Bristol, UK*

<sup>2</sup>*John Innes Centre, Norwich, UK*

<sup>3</sup>*RAGT, Ickleton, Essex, UK*

<sup>4</sup>*Limagrain, Woolpit, Suffolk, UK*

<sup>5</sup>*KWS, Thriplow, Hertfordshire, UK*

<sup>6</sup>*Syngenta Seeds Ltd, Whittlesford, Cambridge, UK*

<sup>7</sup>*KBioscience Unit 7, Hertfordshire, UK*



# Wheat SNP validation – public projects



## Bristol/JIC KASP assays

Cross	Mapped Markers	
Avalon v Cadenza	3,028	1,122 A genome
		1,520 B genome
		386 D genome
Savannah x Rialto	1,543	790 A genome
		594 B genome
		159 D genome
Synthetic x Opata	201	70 A genome
		73 B genome
		49 D genome

- First report of a public linkage map for hexaploid wheat based on KASP to genotype wheat varieties and generate a linkage map.
- 67% polymorphic in varietal screen.
- 4% monomorphic in hexaploid wheat, but polymorphic compared to diploid/tetraploid varieties.

## GCP KASP assays

GCP originator	Validated KASP Assays
Dr Susanne Dreisigacker (CIMMYT)	1,864

- Validation on wheat cultivars originating from Australia, China, India and Mexico.
- Included high anther culture ability and disease resistant varieties.








# Feed the world



# Generation Challenge Program (GCP) panels



# Feed the world



## African Orphan Crops Consortium Genomics Laboratory




- Aim is to reduce stunting by increasing the nutritional value and yields of 100 local African crops by training the best African breeders to use the best tools in the world
- LGC are primary partner of the program and took part in training the 1<sup>st</sup> 25 breeders in Nairobi 2014





# Plant sampling essentials

AOCC Nairobi 2014



# LGC leaf sampling kit



## Kit contents:

- 1 x 96-well tube storage rack with lid
- 12 x perforated strip caps
- 1 x 50g desiccant sachet (in a bag)
- 1 x larger (labelled) sealable bag
- 1 x cutting tool & 1 x cutting mat

If multiple kits have been requested, only 1 x cutting tool and 1 x cutting mat will be sent to you.

We may be able to handle projects outside the specifications given in this document; if this is the case please contact us directly to discuss your requirements:

+44 (0)1992 470757 info.uk@lgcgenomics.com

## Collecting leaf samples

Repeat this process for each plant you want to sample, using a new tube each time.



### Cut leaf discs

- Place the leaf to be sampled on the cutting mat.
- Uncap the cutting tool.
- Cut a disc out of the leaf by pushing the cutting tool into the leaf; twist the tool as you push to make the tool pick the disc up.



### Fill tubes with samples

- Insert the end of the tool into one of the racked tubes, and depress the plunger to dispense the disc.
- *The tubes are in strips of 12, and can be removed from the rack if needed at this step.*
- *Add 1-12 discs per tube, depending on the project (see 'How many disks are required?')*



### Wash cutting tool

- Place the end of the cutting tool in some clean water and depress the plunger a few times.
- Flick the tool until completely dry.

Please ensure the rack is labelled appropriately.

## Preparing samples for transport



### Seal tubes

- Place the perforated Strip caps on top of the tubes.
- Press firmly to ensure caps are secured.
- **Remove the desiccant from its small sealed bag.**



### Add dessicant


- **Place the dessicant sachet on top of the rack of tubes and replace lid.**
- Place completed rack inside the larger (labelled) bag.
- **Force the air out of the bag and then seal it.**



### Prepare for shipping

- Place the sealed bag in a suitable container.
- Our address: **LGC, Extractions department, Units 1 & 2 Trident Industrial Estate, Pindar Road, Hoddesdon, Hertfordshire, UK, EN11 0WZ.**
- *Provide a description of the contents for customs.*

## How many disks are required?





# Agricultural crop improvement case study: The story of the Pigeon pea revolution



- KASP Genotyping assay run by LGC genotyping service lab
- Marker-assisted breeding for improved crops



## Pigeon pea circa. 1992

- Slow growing
- Low yield
- No commercial value
- Neglected by crop breeders
- Un-developed potential



- Breeding program sponsored
- Accelerated MAS using SNPs
- SNPs analysed with KASP

## Pigeon pea circa. 2013

- Drought & disease resistant
- Excellent food supply
- Better market value
- Established life changing market value chain for growers in Tanzania



# Colorado State University Wheat Breeding and Genetics Program



- Scott Healy's group use KASP to identify and monitor traits for Colorado's harsh climate
- They also use KASP to find traits for improving human health such as starch quality and antioxidants.





## *Development of gene specific KASP markers: a toolbox for marker-assisted selection in wheat.*

**Gina Brown-Guedira**<sup>1</sup>, Neelam Kumari<sup>2</sup>, Susan Dreisigacker<sup>3</sup>, Peter Sharp<sup>4</sup>, Catherine Ravel<sup>5</sup>, and Cristobal Uauy<sup>6</sup>.

<sup>1</sup> USDA–ARS Plant Science Research, North Carolina State University, Raleigh, NC, USA; <sup>2</sup> Department of Crop Science, North Carolina State University, Raleigh, NC, USA; <sup>3</sup> International Maize and Wheat Improvement Center (CIMMYT), Mexico; <sup>4</sup> Plant Breeding Institute, University of Sydney, Narellan NSW, Australia; <sup>5</sup> INRA, Université Blaise Pascal, Genetics, Diversity and Ecophysiology of Cereals, Clermont-Ferrand, France; and <sup>6</sup> Department of Crop Genetics, John Innes Centre, Colney, Norwich, United Kingdom.

Whole-genome, SNP detection technologies now available in construction of linkage maps, genome-wide association mapping these technologies are not suited for marker-assisted breeding the same polymorphism on tens of thousands of plants in the same time. A large number of polymorphic markers amenable to high-throughput genotyping are needed and in combination with other marker-assisted breeding approaches. As part of an international collaboration to develop and make publicly available a toolbox of KASP assays for wheat, stream-lined homogeneous assays were developed using a novel assay design, making it well-suited for use in polyploid wheat. Reported sequence variation (SNPs and indels) were used to develop KASP assays for the reduced height, vernalization, and photoperiod-response genes for cloned disease resistance and end-use quality genes. In some cases, the causal gene sequence was known. Also, a number of genes in wheat are associated with important traits. For these, associated SSR and STS markers were evaluated using the Wheat Infinium Assay. SNP in linkage disequilibrium ( $r^2 > 0.9$ ) with the causal gene sequence. In an association mapping approach, markers were developed that are highly predictive for the *Sr36* and *Sbml* resistance genes, and alleles at the *Glu-D1* gene. Association mapping experiments, coupled with SNP genotyping and current sequencing projects in wheat, will result in identification of numerous sequence targets for development of new predictive, homogeneous assays for important genes. We encourage researchers to contribute sequences for development of additional publicly available KASPar assays for use in wheat improvement programs world-wide.

***“KASP provides flexibility in assay design, making it well-suited for use in polyploid wheat, and is low cost”***

***“We encourage researchers to contribute sequences for development of additional publicly available KASPar assays for use in wheat improvement programs world-wide”***

# Plant assay search tool



<b>+</b>	<b>Marker:</b> 3-hydroxy-3-methylglutaryl-coenzyme A reductase 1 <b>Chromosome:</b> 1 <b>Position:</b> 32 <b>Trait:</b> Ggc grey ground colour (seed coat colour)	<a href="#">add to cart</a>
<b>+</b>	<b>Marker:</b> 30S ribosomal protein S31, chloroplastic <b>Chromosome:</b> 1 <b>Position:</b> 32 <b>Trait:</b> Yc cotyledon colour	<a href="#">add to cart</a>
<b>+</b>	<b>Marker:</b> 40S ribosomal protein S15-4 <b>Chromosome:</b> 2 <b>Position:</b> 32 <b>Trait:</b> Ggc grey ground colour (seed coat colour)	<a href="#">add to cart</a>
<b>+</b>	<b>Marker:</b> 40S ribosomal protein S15a-1 <b>Chromosome:</b> 2 <b>Position:</b> 32 <b>Trait:</b> Tgc tan ground colour (seed coat colour)	<a href="#">add to cart</a>
<b>+</b>	<b>Marker:</b> 40S ribosomal protein S18 <b>Chromosome:</b> 2 <b>Position:</b> 32 <b>Trait:</b> Tgc tan ground colour (seed coat colour)	<a href="#">add to cart</a>
<b>+</b>	<b>Marker:</b> 40S ribosomal protein S18 <b>Chromosome:</b> 3 <b>Position:</b> 32 <b>Trait:</b> Tgc tan ground colour (seed coat colour)	<a href="#">add to cart</a>
<b>+</b>	<b>Marker:</b> 50S ribosomal protein L31, chloroplastic <b>Chromosome:</b> 3 <b>Position:</b> 32 <b>Trait:</b> scp seed coat pattern (not close though - 7cM)	<a href="#">add to cart</a>
<b>+</b>	<b>Marker:</b> 60S acidic ribosomal protein P1-2 <b>Chromosome:</b> 3 <b>Position:</b> 32 <b>Trait:</b> Seed diameter, plumpness	<a href="#">add to cart</a>
<b>+</b>	<b>Marker:</b> 60S ribosomal protein L27a-2 <b>Chromosome:</b> 3 <b>Position:</b> 42 <b>Trait:</b> Seed diameter, plumpness	<a href="#">add to cart</a>
<b>+</b>	<b>Marker:</b> ABSCISIC ACID-INSENSITIVE 5-like protein 5 <b>Chromosome:</b> 4 <b>Position:</b> 42 <b>Trait:</b> Seed thickness	<a href="#">add to cart</a>
<b>+</b>	<b>Marker:</b> ABSCISIC ACID-INSENSITIVE 5-like protein 5 <b>Chromosome:</b> 4 <b>Position:</b> 42 <b>Trait:</b> Seed thickness	<a href="#">add to cart</a>

[Reset search](#)

Lentil

Marker Name:

Genetic

Chromosome:

Physical

Mapping Min:

Mapping Max:

Distance/Mapping Interval:

Trait...

**Plant Assay Cart**

[Hide cart](#)

Your cart is empty.

# Web assay search tool

- A web based tool for locating functionally validated SNPs for crops.
- Started out for human validated SNP (more than 100 k!)
- We thought that it might be useful for breeding communities.
- Not there to compete but compliment projects
- We aim for it to:
  - Deposit markers for everyone to use – Can be quickly (minutes) with new markers for sharing.
  - Links to original material
  - Gives basic information about the marker (i.e. sequence or map position (cM))
  - An easy way to order assays

# Assay Search Tool – Submitting

- You can add information linking to reference
- Point me towards where the information is and I'll do the rest.
  - OR
- Provide either map units (cM – reference map) or nucleotide positions (provide the reference genome name).
- Provide marker reference name / sequence
- I'll do the rest



# What will be available?

- Know Pulse – Lentil – core 154 initial assays but more to come
- Sol Cap - All Tomato validated 384 assays and 7k *in silico* designed
- Wheat – cerealsDB assays – around 9000 validated
- More to come!!



# Any Questions?